

CT-WNP

Wrist and Neck Pendant



FEATURES AND BENEFITS

- Specially designed for the Courtney-Thorne 08 Nurse Call System
- Light weight and comfortable so easily accepted by wearers
- Anti-ligature lanyard ensures wearer safety
- Elasticated, adjustable wrist strap for ease of use
- Battery low indicator on both pendant and Nurse Call System ensures early warning of battery replacement
- Wrist strap and lanyard available separately to maximise hygiene
- Location information enabled allowing for efficient attendance
- Spare battery supplied as standard offering a possible 4+ years battery life
- Wipe clean design or full submersion in mild disinfectant for maximum hygiene
- Flashing LED activated on button press for user reassurance
- Reset LED indication giving carer reassurance the unit is reset.
- Configurable call type to manage urgency of response
- Recessed button reduces potential for accidental calls, preventing carers attending residents unnecessarily

SPECIFICATION

Battery life

Over two years operation *

Dimensions

(L) 52.0mm x (W) 44.7mm x (D) 16.0mm

Weight

31g

Call Types

Configurable to Call, Assistance, Emergency, Staff Attack, Reset

Radio frequencies

433MHz or 868MHz

Infra-red frequency

37.9KHz

Enclosure

- Robust polycarbonate plastic with rubberised button
- Elasticated adjustable wrist strap with ergonomic buckle
- Lanyard with anti-ligature protection
- 1 x CR2430 user replaceable coin cell battery
- Spare battery included

Radio transmission range

Up to 100m **

Infra-red transmission range

Up to 10m ***

Standards Compliance

Conforms to relevant EC standards:

- IP67 rated
- EMC
- ROHS/WEEE
- Anti Ligature – Low breaking strain

* *Dependant on features enabled and usage*

** *In free air*

*** *Depending on environment*

Courtney-Thorne

Unit 13
Drewitt Industrial Est.
865 Ringwood Rd
Bournemouth
Dorset
BH11 8LW

Freephone: 0800 068 7419

Email: info@courtney-thorne.co.uk

www.courtney-thorne.co.uk



